
REUSABLE
MANUAL RESUSCITATOR

LATERAL MEDICAL



Flexicare's range of **ergonomically designed** Reusable Manual Resuscitators with **higher levels** of Oxygen delivery provide **fast effective care**, even in the most difficult situations.




Ventilation

Option of 60 or 40cm H_2O Pressure Relief Valves to prevent lung over-inflation and minimise gastric insufflation.

Highest oxygen concentrations at volumes recommended in the ERC 2015 Resuscitation Guidelines1 and International Standards.[2]

SINGLE USE
MANUAL RESUSCITATOR

LATERAL MEDICAL



High oxygen concentrations at optimal tidal volumes, with the option of a blue 60cm H_2O (Adult) or yellow 40cm H_2O (Paediatric & Infant) Pressure Relief Valve, limits excessive ventilation, safeguarding the patient from harm.

Reduced Fatigue
Lightweight construction, low resistance valves and fast recoil bag reduce user fatigue, even during lengthy ventilations.

Grip & Control
Textured surface gives grip and control of the ambidextrous bag, even when wearing gloves and in wet conditions.

Adult Ventilation Bag Strap
Facilitating single handed bag use, the adult bag strap frees the other hand to secure the mask or airway device.

MDI Port
The MDI port is angled so medication is directed into the patient airway, maximising effectiveness and reducing waste.


Optional Manometer
Attaching the single use manometer gives breath-by-breath visualisation of delivered pressures.

Double Swivel Patient End Connector
Valve Assembly rotates 360˚ in two axes for an effective mask seal and bag position, whatever the approach angle.









CapnoCare™

CO_2 MONITORING
ERC 2015 Resuscitation[1] and other Guidelines stipulate the use of capnography during in-hospital or pre-hospital CPR. Monitoring is possible with CapnoCare™ positioned between the Resuscitator and the airway device.

References
1 European Resuscitation Council Guidelines for Resuscitation 2015. Resuscitation, Vol. 95. October 2015.
2 BS EN ISO 10651-4:2009 Lung Ventilators. Part 4: Particular requirements for operator-powered resuscitators.

SINGLE USE
MANUAL RESUSCITATOR

	ADULT			PAEDIATRIC			INFANT			
DESCRIPTION	PART NUMBER	PRESSURE RELIEF VALVE	MASK SIZE	PART NUMBER	PRESSURE RELIEF VALVE	MASK SIZE	PART NUMBER	PRESSURE RELIEF VALVE	MASK SIZE	QTY
Reusable Resuscitator	038-71-950	60cm H_2O	5	038-72-950	40cm H_2O	3	038-73-950	40cm H_2O	1 + 0	5
Reusable Resuscitator + Manometer	038-71-951	60cmH_2O	5	038-72-951	40cm H_2O	3	038-73-951	40cm H_2O	1 + 0	5
Reusable Resuscitator + PEEP Valve	038-71-952	60cm H_2O	5	038-72-952	40cm H_2O	3	038-73-952	40cm H_2O	1 + 0	5
Reusable Resuscitator	038-71-953	-	5	038-72-956	40cm H_2O	2	-	-	-	5
Reusable Resuscitator + Case	038-71-970	60cm H_2O	5	038-72-970	40cm H_2O	3 + 2	038-73-970	40cm H_2O	1 + 0	1

Reusable Replacement Parts

038-71-960	Resuscitator Replacement Patient End Assembly without Pressure Relief Valve	5
038-71-961	Resuscitator Replacement Patient End Assembly with Pressure Relief Valve 60cm H_2O	5
038-71-962	Resuscitator Replacement Patient End Assembly with Pressure Relief Valve 40cm H_2O	5
038-71-965	Resuscitator Intake End Valve Assembly	5
038-71-964	PEEP Valve - Adjustable 5 to 20cm H_2O	5
038-71-967	Resuscitator Replacement Valves, O-Ring & Port Caps Kit	5
038-71-968	Adult Resuscitator Replacement Ventilator Bag with Strap	5
038-72-968	Paediatric Resuscitator Replacement Ventilator Bag	5
038-73-968	Infant Resuscitator Replacement Ventilator Bag	5

Single Use Components

038-71-857	Adult / Paediatric Resuscitator Reservoir Bag with 3m O_2 Tube	10
038-73-857	Infant Resuscitator Reservoir Bag with 3m O_2 Tube	10
038-71-966	Resuscitator Manometer	5
010-750	CapnoCare Capnography Connector with 3m line and Filter with Male Luer	10
010-755	CapnoCare Capnography Connector with 3m line and Filter with Female Luer	10

Single Use Components

038-51-400SIL	Large Adult 5	10
038-51-410SIL	Adult 4	10
038-52-415SIL	Large Child 3	10
038-52-420SIL	Medium Child (Round) 2	10
038-52-430SIL	Small Child / Infant (Round) 1	10
038-53-440SIL	Infant / Neonatal (Round) 0	10

LATERAL MEDICAL

flexicare
TOTAL QUALITY - TOTAL CARE

1300 144 734
sales@lateralmedical.com